

Unlock

Unlock™ is an efficient and gentle cell detachment buffer designed for passaging Axol Neural Progenitor Cells and Axol iPSC-derived Endothelial Cells. It does not contain any mammalian or bacterial components.

Product	Specification
Catalogue number	ax0044
Volume	1 x 25mL
Shipping conditions	Dry Ice
Storage conditions	-80°C
Concentration	1X
Sterility	Sterile
Quality Grade	Cell-culture grade
Xeno-free	Yes

Usage Notes: Thaw overnight (at 4-8°C). Use immediately once thawed. Store for up to 12 months if frozen at - 80°C. Aliquot for use prior to re-freezing.

Advantages

- Free of animal components

Technical Resources

- Application Protocols: <https://www.axolbio.com/page/application-protocols>
- Application Notes: www.axolbio.com/page/application-notes
- Product Page: <https://www.axolbio.com/shop/product/axol-unlock-ax0044-2570?search=ax0044>